

Model 11-T Stretcher



Technical data:

Overall Length: 1915mm

Overall Width: 490mm

Length Folded: 920mm

Depth Folded: 215mm

Height: 180mm

Weight: 10,5kg


Load Capacity: 159kg

Construction materials:

Aluminium Tubing

Aluminium Castings

UV stabilised Flame Resistant Polypropylene



Telephone: +421 327708010 // Fax: +421 327708011

sales@saverman.sk , www.ferno.com

SAVER Manufacturing / Bošáca 893 / 913 07 / Slovakia