




Features

- ARION WHEELED STRETCHER CARRIER is compatible with STANAG 2040 DUO/FOUR NATO STRETCHER and RS6 STRETCHER
- Black sturdy foldable aluminium frame, easy storage
- Can be quickly and easily attached
- 65-0440-001 ARION WHEELED STRETCHER CARRIER

Specifications

Folded Length x Width x Height	Unfolded Length x Width x Height	Wheel diameter	Weight	Load Limit
980 x 640 x 350 mm	1350 x 750 x 880 mm	530 x 50 mm	18 kg	250 kg