



Features

- OXYGEN BAG INTRAXX is bag for a portable oxygen set with the INTRAXX mounting system for attaching the bag in an ambulance or anywhere using the intraxx holder system.
- Made of Complan and Hypalon. The bottom of the bag is made of hypalon, reinforced with abrasion-resistant material. There are 3 plastic hedgehog feet attached to the bottom.
- Contains a single-chamber design that ensures proper segregation of equipment thanks to transparent equipment segregation
- OXYGEN BAG bag has a transport handle for carrying by hand, which is arranged in a way that allows transport in a horizontal position (top center) and the possibility of transport on the shoulder using a detachable, length-adjustable, quick-removable belt
- Numerous reflective elements on the front
- OXYBAG bag is available in two sizes, i.e. size L for transporting cylinders with a capacity of up to 2 L and size XL for transporting cylinders with a capacity of up to 3 L
- The upper part of the bag is equipped with a pocket for marking and identification and an eyelet secured with a flap for leading oxygen and breathing tubes outside.
- The bag's fastening design allows the use of medical equipment (the so-called internal chambers) with the main compartment open, while driving, e.g. in an ambulance, without the risk of its contents falling out.
- Made of abrasion-resistant material, lined with foam ensuring comfort, with the possibility of disinfection, cleaning and washing

60-1800-001 OXYGEN BAG L

60-1800-002 OXYGEN BAG INTRAXX L

60-1800-003 OXYGEN BAG XL

60-1800-004 OXYGEN BAG INTRAXX XL

60-1800-005 OXYGEN BAG XXL + CONNECTOR HOLDER

60-1800-006 OXYGEN BAG XXL INTRAXX + CONNECTOR HOLDER

Specifications

	Height x Width x Depth		Height x Width x Depth	Weight
Size L	30cm x 58cm x 14 cm	Size XL	60cm x 30cm x 17 cm	1,3 kg