

## Model S-240



## Model S-242



### Features

- Designed to be stored in small places
- Folding / Unfolding Chair
- Telescoping front handles with two positions
- Two rear handles (S-240) / Four rear handles (S-242)
- Two rear wheels dia 125 mm (S-240)
- Two wheels dia 80 mm with locks / Two wheels dia 125 mm (S-242)
- Vinyl-coated nylon cover is simple to clean and resists stains, blood and bodily fluids
- Three quick-release patient restraints and leg restraints hold the patient securely

- 60-0110-001 CHAIR S-240
- 60-0110-002 CHAIR S-242

The Models S-240, S-242 Chair has been designed to comply with and exceed international standards such as EN 1865-1:2010 standards.

### Specifications



Length	Width	Height	Weight/without restraints/S-240	Weight/without lower rear handles/S-242	Load limit
74 cm	53 cm	92 cm	8,65 kg	10 kg	170 kg
Length in folded position		Height in folded position		Width in folded position	
92 cm		17 cm		53 cm	