



Features

- Advanced and innovative combat evacuation stretcher
- Meets civilian and military fire safety codes
- Designed for extreme conditions
- Lightweight and sturdy aluminum frame
- Ergonomically designed extendable locking handles
- Slip-resistant mesh, non-flammable
- 6 x IV pole attachment points
- Optional storage bag with backpack straps and carry handles
- Complies with STANAG 2040 and other relevant international standards
- Includes BIO SAFE 4-point restraints
- Meets requirements for roadworthiness (10G-tested)

60-0122-001 STANAG 2040 4-FOLD NATO STRETCHER  NSN 6530400018563

Accessories 65-0317-004 BACKPACK FOR STANAG 2040 4-FOLD STRETCHER

65-0440-001 WHEELED STRETCHER UNDERCARRIAGE ARION

65-1500-050 STANAG 4-FOLD STRETCHER IV POLE



Specifications

Length with extended handles	Length with retracted handles	Width	Height	Folded Length x Width x Height	Weight	Load Capacity
2293 mm	2010 mm	580 mm	145 mm	190 x 230 x 530 mm	8,3 kg	250 kg

Alternatives

60-0122-002 STANAG 2040 FOUR-FOLD NATO STRETCHER with two restraints (does not meet STANAG 2040 Chapter 2 road-worthiness) NSN 6530400018564

60-0102-032 STANAG 2040 DUO-FOLD NATO STRETCHER with two restraints NSN 6530400010645

60-0102-034 STANAG 2040 UNI-FOLD NATO STRETCHER NSN 6530400010644

