



Features

- Independent leg operation provides safe loading into the ambulance, legs unfold automatically for easy unloading
- Color-coded controls are easy to use and require minimal training
- Ambulance loading heights: 500 or 540mm, 560 or 590mm, 650 or 690mm
- Large 200 mm wheels easily roll over any terrain
- Control-end swivel-lock wheels and loading-end fixed wheels with wheel locks
- Four telescoping carrying handles with three positions
- Adjustable backrest with 8 positions up to 90°
- Shock frame and contour frame
- 100% mechanical, low-maintenance cot
- Aluminum construction
- Standard fold-down Sidearms and Foot Support included
- Certifications: EN 1789 and EN 1865

Accessories & Fasteners



Instrument Tray
65-1300-001



IV POLE with extra support
65-1301-001



4- Point Restraint + 2x
Restraint 1pc 180cm
63-4400-001

Mattress org, blk
65-0358-002/90

Mattress welded, blk
65-0358-101

- 60-0303-001 MONOBLOC F2 4P 59CM (H56)
- 60-0323-001 MONOBLOC F2 4P 69CM (H69) FP 3RD POINT SST
- 60-0324-001 MONOBLOC F2 4SW 4P 69CM (H69)
- 60-0304-001 MONOBLOC F2 5P 59CM (H56)
- 60-0325-001 MONOBLOC F2 5P 69CM (H69)



Three-Point Fastener
60-0150-003



Cot Lock F2
60-0253-002

Specifications

Length min	Length max	Width	Height min	Height max	Load Height 1	Load Height 2	Weight	Load Limit	Shipping Weight	Shipping L,W,H
1915 mm	2269 mm	560 mm	272 mm	807 mm	500 mm	540 mm	36 kg	250 kg	43 kg	1910x610x360 mm
					560 mm	590 mm				
					650 mm	690 mm				

