



## Features

- Advanced and innovative combat evacuation stretcher
- Meets civilian and military fire safety codes
- Designed for extreme conditions
- Lightweight and sturdy aluminum frame
- Ergonomically designed extendable locking handles
- Lifeproof vinyl cover with Fire Certification M2/B1, infrared according to Military standards
- Includes 2x restraints
- Optional storage bag with backpack straps and carry handles
- Complies with STANAG 2040 (exception chap. 3.2) and other relevant international standard

60-0102-032 STANAG 2040 2-FOLD NATO STRETCHER      NSN 6530400010645

## Accessories

65-0317-003 BACKPACK NATO STRETCHER

65-0440-001 WHEELED STRETCHER UNDERCARRIAGE



## Specifications

Length with extended handles	Length with retracted handles	Width	Height	Folded Length x Width	Weight	Load Capacity
2290 mm	1988 mm	570 mm	157 mm	937 x 99 mm	8,7 kg	200 kg

## Alternatives

60-0122-001	STANAG 2040 4-FOLD NATO STRETCHER 10G-tested (meets requirements for roadworthiness)	NSN 6530400018563
60-0122-002	STANAG 2040 4-FOLD NATO STRETCHER with two restraints (does not meet STANAG 2040 Chapter 2 road-worthiness)	NSN 6530400018564
60-0102-034	STANAG 2040 UNI-FOLD NATO STRETCHER	NSN 6530400010644

