

## S-265 Scoop Stretcher

Designed so EMS personnel can uncouple either end or both ends of the stretcher and gently scoop up the patient using scissor-type closing motion. Ideal for hip injuries and accident victims found on the street. Supports a patient in the position found, reducing risk of further injury.

Features include a narrow foot-end frame for handling in confined areas, telescoping tubes allow length adjustments to fit patients of various heights, open-center design permits patient to be X-rayed.

- Twin Safety Lock for smooth locking and unlocking
- Blue painted surface (other upon req.)
- Includes three patient restraints

Optional vehicle fastener available.

### S-265 Specifications:

Length: open	166-202 cm
Length: folded	120 cm
Width	42.2 cm
Height: open	6.7 cm
Height: folded	10 cm
Weight	9.75 kg
Loading Capacity	159 kg

